


MOBILE DENTAL ASSISTANT CATALOGUE AND THE PRICE LIST







KLIER
DENTAL

MOBILE DENTAL ASSISTANTS





PRICES EXCLUDING VAT





Assistant 6 drawers		Product code	Standard UNI colors and Wooden decors																								
		MATM6Z																									
		Description	White metal	Antaro	Legrabox	Tandem	Legrabox	Tandem																			
		Drawer depth 400 mm.	1 733 €	1 477 €	1 567 €	1 219 €	1 663 €	1 396 €																			
		Drawers with handles are with soft close. White metal drawers can be filled with dental inserts 3. Antaro, Legrabox and Tandem drawers can be filled with dental inserts 1. Tandem drawers are wooden with part extension runner system.																									
Dimensions		Tandem full extension + 39 € per item.																									
height:	790 mm	Antaro and Legrabox drawers have full extension runner system.																									
depth:	520 mm	The wheels have a brake and a total load capacity of 400 kg.																									
width:	505 mm	Price of the cabinet doesn't include dental inserts.																									
			<table border="1"> <thead> <tr> <th colspan="3">Gloss UNI colors</th> <th colspan="3">Tip-on opening drawers</th> </tr> <tr> <th>Drawers with handles</th> <th>Antaro</th> <th>Legrabox</th> <th>Tandem</th> <th>Legrabox</th> <th>Tandem</th> </tr> </thead> <tbody> <tr> <td>White metal</td> <td>1 906 €</td> <td>1 625 €</td> <td>1 724 €</td> <td>1 340 €</td> <td>1 830 €</td> <td>1 536 €</td> </tr> </tbody> </table>						Gloss UNI colors			Tip-on opening drawers			Drawers with handles	Antaro	Legrabox	Tandem	Legrabox	Tandem	White metal	1 906 €	1 625 €	1 724 €	1 340 €	1 830 €	1 536 €
Gloss UNI colors			Tip-on opening drawers																								
Drawers with handles	Antaro	Legrabox	Tandem	Legrabox	Tandem																						
White metal	1 906 €	1 625 €	1 724 €	1 340 €	1 830 €	1 536 €																					
			<table border="1"> <thead> <tr> <th colspan="3">Matte UNI colors</th> <th colspan="3">Tip-on opening drawers</th> </tr> <tr> <th>Drawers with handles</th> <th>Antaro</th> <th>Legrabox</th> <th>Tandem</th> <th>Legrabox</th> <th>Tandem</th> </tr> </thead> <tbody> <tr> <td>White metal</td> <td>2 079 €</td> <td>1 772 €</td> <td>1 881 €</td> <td>1 462 €</td> <td>1 996 €</td> <td>1 675 €</td> </tr> </tbody> </table>						Matte UNI colors			Tip-on opening drawers			Drawers with handles	Antaro	Legrabox	Tandem	Legrabox	Tandem	White metal	2 079 €	1 772 €	1 881 €	1 462 €	1 996 €	1 675 €
Matte UNI colors			Tip-on opening drawers																								
Drawers with handles	Antaro	Legrabox	Tandem	Legrabox	Tandem																						
White metal	2 079 €	1 772 €	1 881 €	1 462 €	1 996 €	1 675 €																					



Assistant 5 drawers		Product code	Standard UNI colors and Wooden decors																								
		MATM5Z																									
		Description	White metal	Antaro	Legrabox	Tandem	Legrabox	Tandem																			
		Drawer depth 400 mm.	1 693 €	1 438 €	1 528 €	1 179 €	1 624 €	1 357 €																			
		Drawers with handles are with soft close. White metal drawers can be filled with dental inserts 3. Antaro, Legrabox and Tandem drawers can be filled with dental inserts 1. Tandem drawers are wooden with part extension runner system.																									
Dimensions		Tandem full extension + 39 € per item.																									
height:	790 mm	Antaro and Legrabox drawers have full extension runner system.																									
depth:	520 mm	The wheels have a brake and a total load capacity of 400 kg.																									
width:	505 mm	Price of the cabinet doesn't include dental inserts.																									
			<table border="1"> <thead> <tr> <th colspan="3">Gloss UNI colors</th> <th colspan="3">Tip-on opening drawers</th> </tr> <tr> <th>Drawers with handles</th> <th>Antaro</th> <th>Legrabox</th> <th>Tandem</th> <th>Legrabox</th> <th>Tandem</th> </tr> </thead> <tbody> <tr> <td>White metal</td> <td>1 863 €</td> <td>1 581 €</td> <td>1 681 €</td> <td>1 297 €</td> <td>1 786 €</td> <td>1 492 €</td> </tr> </tbody> </table>						Gloss UNI colors			Tip-on opening drawers			Drawers with handles	Antaro	Legrabox	Tandem	Legrabox	Tandem	White metal	1 863 €	1 581 €	1 681 €	1 297 €	1 786 €	1 492 €
Gloss UNI colors			Tip-on opening drawers																								
Drawers with handles	Antaro	Legrabox	Tandem	Legrabox	Tandem																						
White metal	1 863 €	1 581 €	1 681 €	1 297 €	1 786 €	1 492 €																					
			<table border="1"> <thead> <tr> <th colspan="3">Matte UNI colors</th> <th colspan="3">Tip-on opening drawers</th> </tr> <tr> <th>Drawers with handles</th> <th>Antaro</th> <th>Legrabox</th> <th>Tandem</th> <th>Legrabox</th> <th>Tandem</th> </tr> </thead> <tbody> <tr> <td>White metal</td> <td>2 032 €</td> <td>1 725 €</td> <td>1 834 €</td> <td>1 415 €</td> <td>1 949 €</td> <td>1 628 €</td> </tr> </tbody> </table>						Matte UNI colors			Tip-on opening drawers			Drawers with handles	Antaro	Legrabox	Tandem	Legrabox	Tandem	White metal	2 032 €	1 725 €	1 834 €	1 415 €	1 949 €	1 628 €
Matte UNI colors			Tip-on opening drawers																								
Drawers with handles	Antaro	Legrabox	Tandem	Legrabox	Tandem																						
White metal	2 032 €	1 725 €	1 834 €	1 415 €	1 949 €	1 628 €																					

Worktop		Grey with grey core		Black with black core								
White with white core					The worktop is included in the price of the cabinet.							
					Compatible drawer inserts for individual drawer types are marked with the same color.							
					<table border="1"> <tr> <td>Dental inserts for white metal drawers</td> <td></td> </tr> <tr> <td>Dental inserts for the rest types of the drawers</td> <td></td> </tr> </table>						Dental inserts for white metal drawers	
Dental inserts for white metal drawers												
Dental inserts for the rest types of the drawers												
Dimensions		Description										
height:	12 mm	The worktop is made from compact board HPL.										
depth:	505 mm	With groove against the fall of the instruments.										
width:	505 mm											





DENTAL INSERTS 1 - for Antaro, Legrabox, Tandem drawers

ZV1		Product code	ZV2		Product code
		ZV1			ZV2
		Description			Description
		Dental inserts for Antaro, Legrabox and Tandem. Dimensions can be adjusted. (W x D x H): 422 x 457 x 28 mm			Dental inserts for Antaro, Legrabox and Tandem. Dimensions can be adjusted. (W x D x H): 422 x 457 x 28 mm
		PRICE EXCLUDING VAT			PRICE EXCLUDING VAT
		36 €			36 €








ZV3		Product code	ZV4		Product code
		ZV3			ZV4
		Description			Description
		Dental inserts for Antaro, Legrabox and Tandem. Dimensions can be adjusted. (W x D x H): 422 x 457 x 28 mm			Dental inserts for Antaro, Legrabox and Tandem. Dimensions can be adjusted. (W x D x H): 422 x 457 x 57 mm
		PRICE EXCLUDING VAT			PRICE EXCLUDING VAT
		36 €			0 €

ZV5		Product code
		ZV5
		Description
		Dental inserts for Antaro, Legrabox and Tandem. Dimensions can be adjusted. (W x D x H): 422 x 457 x 57 mm
		PRICE EXCLUDING VAT
		36 €

DENTAL ISNERTS 3 - for white metal drawers

ZV740-400		Product code	ZV750-400		Product code
		ZV740-400			ZV750-400
		Description			Description
		Dental insert for white metal drawers 59 mm high. Dimensions cannot be adjusted. (W x D x H): 401 x 382 x 56 mm			Dental insert for white metal drawers 59 mm high. Dimensions cannot be adjusted. (W x D x H): 401 x 382 x 56 mm
		PRICE EXCLUDING VAT			PRICE EXCLUDING VAT
		51 €			57 €

DENTAL ISNERTS 3 - for white metal drawers

<p>ZV751-400</p> 	<p>Product code ZV751-400</p> <p>Description</p> <p>Dental insert for white metal drawers 59 mm high. Dimensions cannot be adjusted. (W x D x H): 401 x 382 x 56 mm</p> <p>PRICE EXCLUDING VAT</p> <p>51 €</p>	<p>ZV752-400</p> 	<p>Product code ZV752-400</p> <p>Description</p> <p>Dental insert for white metal drawers 59 mm high. Dimensions cannot be adjusted. (W x D x H): 401 x 382 x 56 mm</p> <p>PRICE EXCLUDING VAT</p> <p>57 €</p>
<p>ZV753-400</p> 	<p>Product code ZV753-400</p> <p>Description</p> <p>Dental insert for white metal drawers 59 mm high. Dimensions cannot be adjusted. (W x D x H): 401 x 382 x 56 mm</p> <p>PRICE EXCLUDING VAT</p> <p>51 €</p>		
<p>ZV741-400</p> 	<p>Product code ZV741-400</p> <p>Description</p> <p>Dental insert for white metal drawers 120 mm high. Dimensions cannot be adjusted. (W x D x H): 401 x 382 x 117 mm</p> <p>PRICE EXCLUDING VAT</p> <p>59 €</p>	<p>ZV742-400</p> 	<p>Product code ZV742-400</p> <p>Description</p> <p>Dental insert for white metal drawers 120 mm high. Dimensions cannot be adjusted. (W x D x H): 401 x 382 x 117 mm</p> <p>PRICE EXCLUDING VAT</p> <p>59 €</p>
<p>ZV743-400</p> 	<p>Product code ZV743-400</p> <p>Description</p> <p>Dental insert for white metal drawers 120 mm high. Dimensions cannot be adjusted. (W x D x H): 401 x 382 x 117 mm</p> <p>PRICE EXCLUDING VAT</p> <p>57 €</p>	<p>ZV755-400</p> 	<p>Product code ZV755-400</p> <p>Description</p> <p>Dental insert for white metal drawers 120 mm high. Dimensions cannot be adjusted. (W x D x H): 401 x 382 x 117 mm</p> <p>PRICE EXCLUDING VAT</p> <p>51 €</p>

All prices exclude VAT and shipping cost.

Accurate at time of edition 05/2023.

Actual price list available on Klier Dental company website at www.klierdental.eu

KLIER
DENTAL